



## Anti-vibrating tables

- Anti-vibrating tables are built with a sturdy supporting metal base, anti-acid varnished, with adjustable feet and special rubber insulating supports where the anti-vibrating polished marble plate
  - The internal base is protected by an external metal frame, anti-acid varnished, with adjustable feet and complete with buffer guard with front and side panelings avoiding the internal base to be crashed.
  - The front panel is placed backward compared to the working top, allowing in this way the user to work sitting.
  - The working top is lined on all surfaces with HRZ melaminic plastic laminate, white opaque anti-reflex coloured and with rounded edges.
- On request, it is possible to equip the tables with electric facilities series RB 900, with consoles series MB 1000, and lock closure AC 902 for the first upper drawer.

## MESE ANTIVIBRANTE

Cod Produs	Denumirea Produsului		
<u>YP9374</u>	Masa antivibratie cu blat din ceramica 120x75x80H cm cu suport antivibratie . Dim. placa : 48x48 cm		<p>- Suprafata de lucru din folie de plastic mealminata tip HRZ rezistent la acizi si coreziune de culoare opac antireflex, cu orificiu central pentru placa de marmura asezata pe un suport special antivibratie (tip CAU), are deasemenea cadru metalic independent si picioare cu posibilitate de regrare a nivelului</p>
<u>YP9375</u>	Masa antivibratie cu blat din ceramica 120x75x80H cm cu suport antivibratie / <b>cu fundata in nisip</b> Dim. placa : 48x48 cm		
<u>YP9174</u>	Masa antivibratie cu blat din ceramica 60x75x80H cm cu suport antivibratie de Dim. placa : 48x48cm,		
<u>YP9175</u>	Masa antivibratie cu blat din ceramica 60x75x80H cm cu suport antivibratie / <b>cu fundata in nisip</b> Dim. placa : 48x48 cm		
<u>TB 9001</u>	Masa antivibratie cu blat din ceramica 64x75x77H cm cu suport antivibratie Dim. placa : 35x45cm		
<u>TB 9002</u>	Masa antivibratie cu blat din ceramica 120x75x77H cm cu suport antivibratie + dulap cu 4 sertare Dim. placa : 35x45cm		
<u>TB 9003</u>	Masa antivibratie cu blat din ceramica 120x75x77H cm cu suport antivibratie . placa : 35x45cm		
<u>TB 9004</u>	Masa antivibratie cu <b>2 pozitii</b> din blat deceramica 120x75x77H cm cu suport antivibratie Dim. placa : 2 x (35x45cm)		