

Caracteristici tehnice

Performante		Parametre
Partilor Submersibile		Sticla Quartz
Domeniu:	Densitate, g/cm ³ (kg/m ³) Temperatura, °C	0.5000 - 2.0000 (500.0 - 2000.0) + 10 - + 50
Repetabilitate:	densitate, g/cm ³ (kg/m ³) Temperatura, °C	+ 0.0002 + (0.2) + 0.1
Acuratetea	Densitate, g/cm ³ (kg/m ³) Temperatura, °C	± 0.0005 ± (0.5) ± 0.2
Rezolutia	Densitate, g/cm ³ (kg/m ³) Temperatura, °C	± 0.0001 ± (0.1) ± 0.1
Volumul probei		50 ml
Sursa de putere		NiMH baterii 7.2 V 1500 mA/h
Dimensiune (L x W x H)		230 x 130 x 180 mm
Greutate		≈ 2.3 kg
Compensarea Temperaturii		Automat
Afisaj		LCD 2 x 16 caractere