


1




### Creuzet filtrant Filter - crucible

COD CODE	POROZITATE	CAPACITATE	D1
	POROSITY	CAPACITY	
	mm	mL	mm
26.100.1X	0 - 5	8	20
26.100.2X	0 - 5	15	23
26.100.3X	0 - 5	30	30
26.100.4X	0 - 5	50	40

X - va fi inlocuit cu porozitatea solicitata; X - will be replaced with porosity


2



### Adaptor pentru creuzet filtrant Adapter for filter crucible

COD CODE	D1	D2	CREUZET
	mm	mm	mL
26.102.02	27	10	8
26.102.04	34	10	15
26.102.06	42	10	30
26.102.08	50	10	50

3



### Palnie filtrare Filter funnel

COD CODE	CAPACITATE	POROZITATE	D1	D2
	CAPACITY	POROSITY	mm	mm
26.104.1X	50	0 - 5	35	10
26.104.2X	75	0 - 5	45	10
26.104.3X	125	0 - 5	60	10
26.104.4X	500	0 - 5	90	22

X - va fi inlocuit cu porozitatea solicitata; X - will be replaced with porosity


4



### Tub aspiratie Suction tube 1

COD CODE	LUNGIME	DIAMETRUL EXTERIOR
	LENGTH	OUTSIDE DIAMETER
	mm	mm
26.106.02	160	16
26.106.04	180	20

5




### Palnie filtrare Hirsch Suction tube 2

COD CODE	POROZITATE	CAPACITATE	D1	D2	D3
	POROSITY	CAPACITY	mm	mm	mm
26.108.1X	1 - 5	25	55	25	8
26.108.2X	1 - 5	70	100	30	10

X - va fi inlocuit cu porozitatea solicitata; X - will be replaced with porosity

6



### Palnie conducta Pipeline filter

COD CODE	POROZITATE	D1	D2
	POROSITY	mm	mm
26.110.1X	0 - 5	12	30
26.110.2X	0 - 5	14	60
26.110.3X	0 - 5	18	90

X - va fi inlocuit cu porozitatea solicitata; X - will be replaced with porosity





