



TYPE 34S

Large Cylindrical Polarimeter Cell with Screw Caps

- Includes 2 SC2 open screw caps and two SC1 closed screw caps

L/P (mm)	Outer Dimensions (mm)			Inside Dimensions (mm)		Material & Part Numbers	
	LENGTH	HEIGHT	DIAMETER	DIAMETER	LENGTH	CAPACITY (ml)	UV
50	52.5	35.5	22.0	19.0	50.0	14.000	34SUV50
100	102.5	35.5	22.0	19.0	100.0	28.000	34SUV100




TYPE 41

Standard cuvette with screw cap

- Two windows polished
- Height includes threading, but not screw cap
- Cuvette comes with both open and closed screw caps
- Each cap contains a PTFE septa seal
- Suitable for anaerobic work

L/P (mm)	Outer Dimensions (mm)			Inside Dimensions (mm)		Material & Part Numbers		
	LENGTH	WIDTH	HEIGHT	LENGTH	WIDTH	CAPACITY (ml)	G	UV
1	3.5	12.5	55.0	1.0	10.0	0.400	-	41UV1
2	4.5	12.5	55.0	2.0	10.0	0.700	-	41UV2
5	7.5	12.5	55.0	5.0	10.0	1.700	41G5	41UV5
10	12.5	12.5	55.0	10.0	10.0	3.500	41G10	41UV10
50	52.5	12.5	55.0	50.0	10.0	17.500	41G50	41UV50
100	102.5	12.5	55.0	100.0	10.0	35.000	41G100	41UV100

SCREW CAP OD: 16mm
SEPTA: Diameter: 11mm
Thickness: 3.5mm
THREADING ID: 8mm
THREADING OD: 12mm



TYPE 46

Variable capacity micro/semi-micro cuvette with screw cap

- Two windows polished
- Height includes threading, but not screw cap
- Cuvette comes with both open and closed screw caps
- Each cap contains a PTFE septa seal
- Suitable for anaerobic work
- Base thickness is 3mm

L/P (mm)	Outer Dimensions (mm)			Inside Dimensions (mm)		Material & Part Numbers	
	LENGTH	WIDTH	HEIGHT	LENGTH	WIDTH	CAPACITY (ml)	UV
10	12.5	12.5	55.0	10.0	1.0	0.350	46UV10X1
10	12.5	12.5	55.0	10.0	2.0	0.700	46UV10X2
10	12.5	12.5	55.0	10.0	3.0	1.050	46UV10X3
10	12.5	12.5	55.0	10.0	4.0	1.400	46UV10X4
10	12.5	12.5	55.0	10.0	5.0	1.750	46UV10X5

SCREW CAP OD: 16mm
SEPTA: Diameter: 11mm
Thickness: 3.5mm
THREADING ID: 8mm
THREADING OD: 12mm