



MultiLab

B-190 Series - Accessories

ACCESSORIES FOR B-190 SERIES

M-001	Eyepiece H5x.
M-002.1	Eyepiece WF10x/18mm.
M-003	Eyepiece WF16x/12mm.
M-004	Eyepiece with micrometer WF10x/18mm, 10mm/100um.
M-008	Eyepiece with pointer WF10x/18mm.
M-137	ACH (Achromatic) objective 4x/0.10
M-138	ACH (Achromatic) objective 10x/0.25.
M-139	ACH (Achromatic) objective 20x/0.40.
M-141	ACH (Achromatic) objective 40x/0.65.
M-142	ACH (Achromatic) objective 60x/0.80.
M-143	ACH (Achromatic) objective 100x/1.25 (Oil).
M-190	Polarizing set (filters only).
DC-002	Plastic dust cover, medium.
M-069	Solar battery pack.
M-974	Blue filter, 32 mm diameter.
M-976	Green filter, 32 mm diameter.
M-978	Yellow filter, 32 mm diameter.
M-988	Frosted glass filter, 32 mm diameter.
M-971	Plane-concave mirror, with base.
M-114	C-Mount adapter for 1/2" sensor.
M-115	C-Mount adapter for 1/3" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
15008	OPTIKA immersion oil 10ml.
M-005	Micrometer slide for software calibration, 1mm/10um, 10mm/100um.
15104	Cleaning kit.
VP-190	IQ/OQ/PQ Validation Protocols.

M-069 - Solar battery pack

Included battery: rechargeable – Lithium-Poly. Capacity: 2600 mAh.
Output voltage: 5,5Vdc. - Dimensions: 120x73x10 mm.
Autonomy: over 6 hours at medium intensity (X-LED³).
Charging modes: with solar panel (12h),
with external USB power supply (not included)
or from PC USB port (5h).



B-192 Powered by solar battery pack (M-069)

15104 - Cleaning kit

It cleans glass quickly and effectively,
without leaving residue or odor.
Ideal for precision lens or prism cleaning.



Multi Lab SRL

Str. Baia de Arama, nr.1, Bloc C
Spatiu Comercial Nr. 2,
Sector 2 - Cod.Postal 022204
Bucuresti - Romania

Tel.+4 021 3127144 / 45 :: Fax.+4 021 3127141

www.multilab.ro :: www.labshop.ro :: www.eulab.eu