

B-290 Series - Accessories

ACCESSORIES FOR B-290 SERIES

M-160	Eyepiece WF10x/20mm.
M-161	Eyepiece WF15x/16mm.
M-162	Eyepiece WF20x/10mm.
M-163	Micrometer eyepiece WF10x/20mm, 10mm/100um.
M-008.1	Eyepiece with pointer WF10x/20mm.
M-164	E-PLAN objective 4x/0.10 (only for B-292 / B-293 / B-290TB).
M-165	E-PLAN objective 10x/0.25 (only for B-292 / B-293 / B-290TB).
M-166	E-PLAN objective 20x/0.40 (only for B-292 / B-293 / B-290TB).
M-167	E-PLAN objective 40x/0.65 (only for B-292 / B-293 / B-290TB).
M-168	E-PLAN objective 60x/0.80 (only for B-292 / B-293 / B-290TB).
M-169	E-PLAN objective 100x/1.25 (Oil) (only for B-292 / B-293 / B-290TB).
M-144	IOS E-PLAN objective 4x/0.10 (only for B-292PLi / B-293PLi).
M-145	IOS E-PLAN objective 10x/0.25 (only for B-292PLi / B-293PLi).
M-146	IOS E-PLAN objective 20x/0.40 (only for B-292PLi / B-293PLi).
M-147	IOS E-PLAN objective 40x/0.65 (only for B-292PLi / B-293PLi).
M-148	IOS E-PLAN objective 100x/1.25 (Oil) (only for B-292PLi / B-293PLi).
M-149	IOS E-PLAN objective 60x/0.80 (only for B-292PLi / B-293PLi).
M-174	Polarizing set, filters only (for B-290 series).
M-184	Darkfield condenser.
M-175	Rotating table for polarizing set.
M-069	Solar battery pack.
M-971	Plane-concave mirror, with base.
M-114	C-Mount adapter for 1/2" sensor.
M-115	C-Mount adapter for 1/3" sensor.
M-116	C-Mount adapter for 2/3" sensor.
M-173	Photo adapter for APS-C and Full Frame Reflex cameras.
M-975	Blue filter, 45mm diameter.
M-977	Green filter, 45mm diameter.
M-979	Yellow filter, 45mm diameter.
M-989	Frosted glass filter, 45mm diameter.
M-666.290	Heating stage, with digital temperature controller.
M-005	Micrometer slide for software calibration, 1mm/10um, 10mm/100um.
M-031.1	Plastic dust cover – For 290 Series.
M-031.3	TNT dust cover – For B-290TB.
15008	OPTIKA immersion oil 10ml.
15104	Cleaning kit.
VP-290	IQ/OQ/PQ Validation Protocols.

M-069 - Solar battery pack

Included battery: rechargeable – Lithium-Poly. Capacity: 2600 mAh.
 Output voltage: 5.5Vdc. - Dimensions: 120x73x10 mm.
 Autonomy: over 6 hours at medium intensity (X-LED³).
 Charging modes: with solar panel (12h),
 with external USB power supply (not included)
 or from PC USB port (5h).



15104 - Cleaning kit

It cleans glass quickly and effectively,
 without leaving residue or odor.
 Ideal for precision lens or prism cleaning.

