

B-380 Series - Accessories

ACCESSORIES FOR B-380

- M-160 Eyepiece WF10x/20mm.
- M-161 Eyepiece EW15x/16mm.
- M-162 Eyepiece WF20x/10mm.
- M-163 Micrometer eyepiece WF10x/20mm, 10mm/100um.
- M-008.1 Eyepiece with pointer WF10x/20mm.
- M-164 E-PLAN objective 4x/0.10 (only for B-382PL-ALC, B-383PL, B-383DK).
- M-165 E-PLAN objective 10x/0.25 (only for B-382PL-ALC, B-383PL, B-383DK).
- M-166 E-PLAN objective 20x/0.40 (only for B-382PL-ALC, B-383PL, B-383DK).
- M-167 E-PLAN objective 40x/0.65 (only for B-382PL-ALC, B-383PL, B-383DK).
- M-168 E-PLAN objective 60x/0.80 (only for B-382PL-ALC, B-383PL, B-383DK).
- M-169 E-PLAN objective 100x/1.25 (Oil) (only for B-382PL-ALC, B-383PL).
- M-144 IOS E-PLAN objective 4x/0.10 (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL).
- M-145 IOS E-PLAN objective 10x/0.25 (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL).
- M-146 IOS E-PLAN objective 20x/0.40 (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL).
- M-147 IOS E-PLAN objective 40x/0.65 (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL).
- M-149 IOS E-PLAN objective 60x/0.80 (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL).
- M-148 IOS E-PLAN objective 100x/1.25 (Oil) (Oil) (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL).
- M-059 PLAN objective 100x/1.25OI (oil iris) (only for B-383DK).
- M-170 PLAN PH objective 10x/0.25 (only for B-383PH, B-382PH-ALC).
- M-171 PLAN PH objective 20x/0.40 (only for B-383PH, B-382PH-ALC).
- M-172 PLAN PH objective 40x/0.65 (only for B-383PH, B-382PH-ALC).
- M-182 PLAN PH objective 100x/1.25 (Oil) (only for B-383PH, B-382PH-ALC).
- M-183 IOS PLAN PH objective 10x/0.25 (only for B-383PHi, B-382PHi-ALC).
- M-176 IOS PLAN PH objective 20x/0.40 (only for B-383PHi, B-382PHi-ALC).
- M-177 IOS PLAN PH objective 40x/0.65 (only for B-383PHi, B-382PHi-ALC).
- M-178 IOS PLAN PH objective 100x/1.25 (Oil) (only for B-383PHi, B-382PHi-ALC).
- M-144P IOS E-PLAN POL objective 4x/0.10 (only for B-383POL).
- M-145P IOS E-PLAN POL objective 10x/0.25 (only for B-383POL).
- M-146P IOS E-PLAN POL objective 20x/0.40 (only for B-383POL).
- M-147P IOS E-PLAN POL objective 40x/0.65 (only for B-383POL).
- M-149P IOS E-PLAN POL objective 60x/0.80 (only for B-383POL).
- M-148P IOS E-PLAN POL objective 100x/1.25 (Oil) (only for B-383POL).
- M-335 IOS PLAN MET objective 50x/0.75 (only for B-383LD1, B-383LD2, B-383FL, B-383MET).
- M-337 IOS PLAN MET objective 4x/0.10 (only for B-383MET).
- M-338 IOS PLAN MET objective 10x/0.25 (only for B-383MET).
- M-339 IOS PLAN MET objective 20x/0.40 (only for B-383MET).
- M-179 Complete PH Set – PLAN PH obj. 10x, 40x, 100x - Condenser for BF/DF/PH (only for B-383PL, B-382PL-ALC, B-383DK).
- M-181 Complete PH Set – IOS PLAN PH obj. 10x, 20x, 40x, 100x – Condenser for BF/DF/PH (only for B-382PLi-ALC, B-383PLi, B-383LD1, B-383LD2, B-383FL, B-383MET).
- M-151 HBO100W high-pressure mercury bulb for fluorescence (for B-383FL only).
- M-666.380 Heating stage, with digital temperature controller.
- M-174.1 Polarizing set, filters only (for B-380 series).
- M-175 Rotating table for polarizing set.
- M-185 DF condenser for dry objectives.
- M-975.1 Ring with blue filter, 45mm diameter.
- M-977.1 Ring with green filter, 45mm diameter.
- M-979.1 Ring with yellow filter, 45mm diameter.
- M-989.1 Ring with frosted glass filter, 45mm diameter.
- M-173 Photo adapter for APS-C and Full Frame Reflex cameras.
- M-114 C-Mount adapter for 1/2" sensor.
- M-115 C-Mount adapter for 1/3" sensor.
- M-116 C-Mount adapter for 2/3" sensor.
- M-005 Micrometer slide for software calibration, 1mm/10um, 10mm/100um.
- DC-003 TNT dust cover, medium – 600(l)x550(h)mm.
- 15008 OPTIKA immersion oil 10ml.
- 15104 Cleaning kit.
- VP-380 IQ/OQ/PQ Validation Protocols.

M-069 - Solar battery pack

Included battery: rechargeable – Lithium-Poly. Capacity: 2600 mAh.
Output voltage: 5.5Vdc. - Dimensions: 120x73x10 mm.
Autonomy: over 6 hours at medium intensity (X-LED³).
Charging modes: with solar panel (12h),
with external USB power supply (not included)
or from PC USB port (5h).



15104 - Cleaning kit

It cleans glass quickly and effectively,
without leaving residue or odor.
Ideal for precision lens or prism cleaning.



MultiLab SRL
Tel. +4 021 312 7144 / 45
contact@multilab.ro :: www.multilab.ro